


MERCODIA OXIDIZED LDL ELISA

10-1143-01



WHEN ACCURACY MATTERS

www.mercodia.com

Intended use

A high quality enzyme immunoassay for the quantification of human oxidized LDL in serum or plasma.

Test principle

Mercodia Oxidized LDL ELISA is a solid phase two-site enzyme immunoassay based on the sandwich technique, in which two monoclonal antibodies are directed against separate antigenic determinants on the oxidized apolipoprotein B molecule. Oxidized LDL in the sample reacts with anti-oxidized LDL antibodies bound to microtitration wells and peroxidase-conjugated anti-oxidized LDL antibodies in the solution.

Summary of protocol

- Add 25 µL Calibrators, controls and samples
- Add 100 µL Assay Buffer
- Incubate 2 hours on shaker at room temperature
- Wash plate 6 times
- Add 100 µL enzyme conjugate solution
- Incubate 1 hour on shaker at room temperature
- Wash plate 6 times
- Add 200 µL Substrate TMB
- Incubate 15 minutes at room temperature
- Add 50 µL Stop Solution
- Shake for approximately 5 seconds on shaker
- Measure at 450 nm

Samples

The recommended use of specimen in the Mercodia Oxidized LDL ELISA is fresh EDTA-plasma. Serum or heparin plasma may also be used. Samples are diluted 1/6561 in a two-step procedure in Sample Buffer:

1. Sample 25 µL
Sample Buffer 2000 µL
Cap tubes, invert 3 times and vortex Dilution 1/81
2. 1/81 dilution 25 µL
Sample buffer 2000 µL
Cap tubes, invert 3 times and vortex Dilution 1/6561

Catalog No

10-1143-01 1 x 96 wells

Measurement range

1.4 - 22.5 mU/L (8 - 150 U/l when multiplied with dilution factor)

Test characteristics

Sensitivity

The detection limit is < 1 mU/L calculated as three standard deviations above the Calibrator 0.

Recovery

Recovery upon addition is 85 - 107 % (95 %).

Precision

Each sample was analyzed in 3 - 8 replicates on 20 different occasions.

Sample	Mean value (mU/L)	Coefficient of variation		
		within assay %	between assay %	total assay %
1	8.5	5.5	6.2	8.3
2	19	7.3	4.0	8.3
3	32	6.2	4.0	7.4

Performance limitations

Grossly lipemic, icteric or haemolyzed samples do not interfere in the assay.

Specificity

Specific for human oxidized LDL.

Global Headquarters

Mercodia AB
Phone +46 18 57 00 70
E-mail info@mercodia.com

US Sales Office

Mercodia Inc
Phone (336) 725-8623
E-mail info_usa@mercodia.com